

Thermo Scientific Dosimetry Solution provides a simple and convenient solution to protect your most valuable assets, people and facilities, while recording and tracking radiation dose measurements for peace of mind.

Thermo Scientific Dosimetry Solution

Dosimetry Solution Kit




Key Features

- One kit provides total dosimetry solution
- Advanced radiological performance
- Proven technology
- Small and very lightweight dosimeter
- Low power consumption, EPD operates on a single battery
- Rugged ergonomic design
- Integrated infrared reader in carry case
- EasyEPD and Easy Issue software included

Applications

- Medical walk-in clinics
- Radiology labs and stand-alone facilities
- Veterinarian offices
- Dental offices
- Ambulance companies
- First responders
- Research laboratories



The Thermo Scientific Dosimetry Solutions kit is designed for providing a complete and comprehensive dosimetry solution to facilities that do not need or cannot afford a large dosimetry system. Building on the proven technology of the Thermo Scientific Mk2 Electronic Personal Dosimeter (EPD), the standard kit contains 5 EPDs, software to configure and issue the EPDs, and a rugged carrying case with an integrated IrDA reader.

The convenience of the Dosimetry Solutions kit allows administrators of clinics, medical offices, labs, and other locations to receive all the dosimetry equipment needed in one package. There is no need to configure a system from a variety of components. Instead, the Dosimetry Solution reduces any possible confusion surrounding how to accurately measure radiation dose rates.

Depending on the need of the customer, a kit can contain dosimeters to measure for beta/gamma, neutron or gamma only. For larger facilities, there is an option to order from 5 to up to 12 Electronic Personal Dosimeters. No experience in dosimetry is required to operate the dosimeters, meaning you save costs on training.

The high quality of the Thermo Scientific Mk2 EPD, combined with the included Easy Issue and EasyEPD2 software packages and the integrated infrared EPD reader, makes the Dosimetry Solution the perfect kit for any facility looking to record dose rates. This provides the facility with insurance against possible future litigation or regulatory actions.

Components of the Dosimetry Solution kit

Electronic Personal Dosimeter

- Mk2, N or G dosimeter models can be selected; standard kit contains 5 EPDs, option to select up to 12 EPDs
- Multi-detector technology
- Excellent response to all radiation types; choice of beta/gamma or gamma-only EPD
- Lightweight design: 95 g (3.2 oz) including battery and clip
- Size: 85 mm x 63 mm x 19 mm (3.3" x 2.5" x .08"), excluding clip
- Power supply: single AA battery 1.5 V alkaline cell for up to 30 weeks of continuous operation OR 3.6 V lithium battery for typically 20-40 weeks continuous operation (interchangeable; life estimates depend on type)
- 10-year data retention without battery
- Configurable audible and visual alarms for dose, dose rate, count down time, read time, and failure modes
- Excellent immunity to electromagnetic interference
- Easy-to-read display with optional backlighting

EasyEPD2 Software

- No database required
- Configuration software to customize alarm levels, display information, and more
- Interfaces directly to the Thermo IrDA reader for the Mk2 EPD
- Automatic "discovery" and reading of a detected EPD unit
- Basic EPD Issue and Return functions supported
- Individual screen windows for each logical functional area of control
- Compatible with Windows™

Easy Issue Software

- Flexible, simple and scalable database that can be stand-alone or networkable
- Supports grab-and-go operations
- Ability to allocate alarm levels for different tasks
- Restrict EPD issue based on personal radiation dose levels

EPD IrDA Reader

- Infrared EPD reader integrated into kit carrying case
- Reader software included
- Connect directly to a PC via serial interface
- Requires no power supply

Carrying Case

- Rugged Pelican case with molded foam inserts to protect EPD units
- Integrated EPD IrDA reader in case
- Easy storage for the entire system in one self-contained unit
- Can be shipped using FedEx or UPS

©2009 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Results may vary under different operating conditions. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representatives for details. Literature Code RMP DIAB 2009

Worldwide
Frauenauracher Strasse 96 +49 (0) 9131 909-0
D 91056 Erlangen, Germany +49 (0) 9131 909-205 fax

United Kingdom
Bath Road, Beenham, +44 (0) 118 971 2121
Reading RG7 5PR United Kingdom +44 (0) 118 971 2835 fax

United States +1 (508) 520-2815
27 Forge Parkway +1 (800) 274-4212 toll-free
Franklin, MA 02038 USA +1 (508) 428-3535 fax

www.thermo.com/rmp

Thermo
SCIENTIFIC